

# Package: racir (via r-universe)

October 31, 2024

**Type** Package

**Title** Rapid A/Ci Response (RACiR) Data Analysis

**Version** 2.0.0

**Date** 2020-07-29

**Maintainer** Joseph Stinziano <josephstinziano@gmail.com>

**Description** Contains functions useful for reading in Licor 6800 files, correcting and analyzing rapid A/Ci response (RACiR) data. Requires some user interaction to adjust the calibration (empty chamber) data file to a useable range. Calibration uses a 1st to 5th order polynomial as suggested in Stinziano et al. (2017) <doi:10.1111/pce.12911>. Data can be processed individually or batch processed for all files paired with a given calibration file. RACiR is a trademark of LI-COR Biosciences, and used with permission.

**URL** <https://github.com/jstinzi/racir>

**BugReports** <https://github.com/jstinzi/racir/issues>

**License** MIT + file LICENSE

**Depends** R (>= 4.0.0)

**Encoding** UTF-8

**LazyData** true

**RoxygenNote** 7.1.1

**Suggests** knitr, rmarkdown, testthat

**VignetteBuilder** knitr

**Imports** utils, stats, graphics

**Repository** <https://jstinzi.r-universe.dev>

**RemoteUrl** <https://github.com/jstinzi/racir>

**RemoteRef** HEAD

**RemoteSha** 3850ebf0fa0e475c25b3798fbb8fe7da1fb37839

## Contents

racircal . . . . .	2
racircalbatch . . . . .	3
racircalbatch_advanced . . . . .	4
racircalcheck . . . . .	5
racircalcheck_advanced . . . . .	6
racircal_advanced . . . . .	7
read_6800 . . . . .	9

<b>Index</b>	<b>10</b>
--------------	-----------

---

racircal	<i>Corrects rapid A/Ci response (RACiR) data from leaves using empty chamber data.</i>
----------	--

---

### Description

racircal Corrects your RACiR data based on calibration data. Produces corrected A vs. Ci graph. Output is a data frame with corrected RACiR data using variable names Acor and Cicor for the corrected A and Ci values.

### Usage

```
racircal(
  data,
  caldata,
  mincut,
  maxcut,
  title,
  varnames = list(A = "A", Ca = "Ca", CO2_r = "CO2_r", E = "E", gtc = "gtc")
)
```

### Arguments

data	Data frame with the RACiR response data
caldata	Data frame with the calibration data
mincut	Minimum cutoff value for reference CO2 (CO2_r). Used to cut out the data from the initial chamber mixing. Default value is set to the minimum COR_r value.
maxcut	Maximum cutoff value for reference CO2 (CO2_r). Used to cut out the data from the end of the response. Not needed in all cases. Default value is set to the maximum COR_r value.
title	Title of output graph - useful for batch RACiR corrections.
varnames	Variable names - this allows for the use of this code with other machines and setups where variable names may differ.

**Value**

racircal returns a data frame with corrected RACiR data

**Examples**

```
#Read in data
data <- read_6800(system.file("extdata", "poplar_2", package = "racir"))
caldata <- read_6800(system.file("extdata", "cal", package = "racir"))
#Correct data
data_corrected <- racircal(data = data, caldata = caldata,
                           mincut = 350, maxcut = 780, title = "Test")
```

---

racircalbatch	<i>Corrects a batch of rapid A/Ci response (RACiR) data</i>
---------------	---

---

**Description**

racircalbatch Corrects your RACiR data files based on a calibration file. Produces diagnostic graphs of A vs. Ci for quality control. Output includes a list of data frames with corrected data.

**Usage**

```
racircalbatch(
  caldata,
  data,
  mincut,
  maxcut,
  title,
  varnames = list(A = "A", Ca = "Ca", CO2_r = "CO2_r", E = "E", gtc = "gtc")
)
```

**Arguments**

caldata	Data frame with the calibration data
data	List of data frames with the RACiR response data
mincut	Minimum cutoff value for reference CO <sub>2</sub> (CO <sub>2_r</sub> ). Used to cut out the data from the initial chamber mixing. Default value is set to the minimum COR <sub>r</sub> value.
maxcut	Maximum cutoff value for reference CO <sub>2</sub> (CO <sub>2_r</sub> ). Used to cut out the data from the end of the response. Not needed in all cases. Default value is set to the maximum COR <sub>r</sub> value.
title	Vector for titles of output graph - useful for batch RACiR corrections. Length must be equal to data list length
varnames	Variable names - this allows for the use of this code with other machines and setups where variable names may differ.

**Value**

racircalbatch calibrates a batch of RACiR data

**Examples**

```
#Create a list of files
files <- c(system.file("extdata", "poplar_1", package = "racir"),
           system.file("extdata", "poplar_2", package = "racir"))
data <- vector("list", length(files))
for(i in seq_along(files)){
  data[[i]] <- read_6800(files[i])
  names(data)[i] <- files[i]
}

caldata <- read_6800(system.file("extdata", "cal", package = "racir"))
output <- racircalbatch(caldata = caldata, data = data,
                       mincut = 300, maxcut = 780, title = files)
```

---

racircalbatch\_advanced

*Corrects a batch of rapid A/Ci response (RACiR) data*

---

**Description**

racircalbatch\_advanced Corrects your RACiR data files based on a calibration file. Produces diagnostic graphs of A vs. Ci for quality control. Output includes a list of data frames with corrected data.

**Usage**

```
racircalbatch_advanced(
  caldata,
  data,
  mincut,
  maxcut,
  digits,
  title,
  varnames = list(A = "A", Ca = "Ca", CO2_r = "CO2_r", E = "E", gtc = "gtc")
)
```

**Arguments**

caldata	Data frame with the calibration data
data	List of data frames with the RACiR response data

mincut	Minimum cutoff value for reference CO2 (CO2_r). Used to cut out the data from the initial chamber mixing. Default value is set to the minimum COR_r value.
maxcut	Maximum cutoff value for reference CO2 (CO2_r). Used to cut out the data from the end of the response. Not needed in all cases. Default value is set to the maximum COR_r value.
digits	Specifies rounding for groups. Defaults to -2 (100s). Effectively uses 100 ppm intervals (e.g. data matching >50 ppm to 150 ppm would be assigned to an interval centered around 100 ppm for reference CO2).
title	Vector for titles of output graph - useful for batch RACiR corrections. Length must be equal to data list length
varnames	Variable names - this allows for the use of this code with other machines and setups where variable names may differ.

### Value

racircalbatch\_advanced uses racircal\_advanced on many files

### Examples

```
#Create a list of files
files <- c(system.file("extdata", "poplar_1", package = "racir"),
           system.file("extdata", "poplar_2", package = "racir"))
data <- vector("list", length(files))
for(i in seq_along(files)){
  data[[i]] <- read_6800(files[i])
  names(data)[i] <- files[i]
}

caldata <- read_6800(system.file("extdata", "cal", package = "racir"))
output <- racircalbatch_advanced(caldata = caldata, data = data,
                                mincut = 300, maxcut = 780, title = files)
```

---

racircalcheck	<i>Allows visual checking of rapid A/Ci response (RACiR) calibration data using empty chamber data.</i>
---------------	---

---

### Description

racircalcheck Used to check range of calibration file. Produces diagnostic graphs of A vs. Ci for quality control. Output includes plots for checking and confirming cutoff values, and a plot with the fit, as well as information as to which polynomial fit the data best.

**Usage**

```
racircalcheck(
  data,
  mincut,
  maxcut,
  varnames = list(A = "A", Ca = "Ca", CO2_r = "CO2_r", E = "E", gtc = "gtc")
)
```

**Arguments**

data	Data frame with the calibration (empty chamber) rapid A/Ci response
mincut	Minimum cutoff value for reference CO <sub>2</sub> (CO <sub>2_r</sub> ). Used to cut out the data from the initial chamber mixing. Default value is set to the minimum COR_r value.
maxcut	Maximum cutoff value for reference CO <sub>2</sub> (CO <sub>2_r</sub> ). Used to cut out the data from the end of the response. Not needed in all cases. Default value is set to the maximum COR_r value.
varnames	Variable names - this allows for the use of this code with other machines and setups where variable names may differ.

**Value**

racircalcheck allows visual checking of RACiR calibration data

**Examples**

```
#Read in the file
data <- read_6800(system.file("extdata", "cal", package = "racir"))
#Run calibration check
racircalcheck(data = data,
              mincut = 350,
              maxcut = 780)
```

---

racircalcheck\_advanced

*Allows visual checking of rapid A/Ci response (RACiR) calibration data using empty chamber data.*

---

**Description**

racircalcheck\_advanced Used to check range of calibration file. Produces diagnostic graphs of A vs. Ci for quality control.

**Usage**

```
racircalcheck_advanced(
  data,
  mincut,
  maxcut,
  digits,
  varnames = list(A = "A", Ca = "Ca", CO2_r = "CO2_r", E = "E", gtc = "gtc")
)
```

**Arguments**

data	Data frame with the calibration (empty chamber) rapid A/Ci response
mincut	Minimum cutoff value for reference CO <sub>2</sub> (CO <sub>2_r</sub> ). Used to cut out the data from the initial chamber mixing. Default value is set to the minimum COR_r value.
maxcut	Maximum cutoff value for reference CO <sub>2</sub> (CO <sub>2_r</sub> ). Used to cut out the data from the end of the response. Not needed in all cases. Default value is set to the maximum COR_r value.
digits	Specifies rounding for groups. Defaults to -2 (100s). Effectively uses 100 ppm intervals (e.g. data matching >50 ppm to 150 ppm would be assigned to an interval centered around 100 ppm for reference CO <sub>2</sub> ).
varnames	Variable names - this allows for the use of this code with other machines and setups where variable names may differ.

**Value**

racircalcheck\_advanced returns a data frame with corrected RACiR data

**Examples**

```
#Read in data
data <- read_6800(system.file("extdata", "poplar_2", package = "racir"))
caldata <- read_6800(system.file("extdata", "cal", package = "racir"))
#Correct data
racircalcheck_advanced(data = data, mincut = 350, maxcut = 780)
```

---

racircal_advanced	<i>Corrects rapid A/Ci response (RACiR) data from leaves using empty chamber data.</i>
-------------------	--

---

**Description**

racircal\_advanced Interval correction for RACiR data.





---

read_6800	<i>Reads files from the Li-Cor 6800</i>
-----------	---

---

**Description**

read\_6800 Reads Li-Cor 6800 files, which are delimited by spaces and tabs.

**Usage**

read\_6800(x)

**Arguments**

x                    A Li-Cor 6800 data file name of the form: "mydata".

**Value**

read\_6800 imports a Li-Cor 6800 file as a data frame

# Index

racircal, [2](#)  
racircal\_advanced, [7](#)  
racircalbatch, [3](#)  
racircalbatch\_advanced, [4](#)  
racircalcheck, [5](#)  
racircalcheck\_advanced, [6](#)  
read\_6800, [9](#)